



# CPX™ 4 Breast Tissue Expander Injection Port

The dimensions below are measurements of the injection port for the CPX™ 4 Tissue Expander.



## CPX™ 4 SILTEX™ Breast Tissue Expander

**Injection Port**  
see detailed  
illustration at top



**MRI Unsafe.** The injection port contains a stainless steel magnetic port, which may be identified during a metal detection procedure.

### Important Safety Information:

MENTOR® CPX™ 4 Breast Tissue Expanders can be utilized for breast reconstruction after mastectomy, correction of an underdeveloped breast, scar revision and tissue defect procedures. These expanders are intended for temporary subcutaneous or submuscular implantation; they should be used within a time frame determined by the physician to achieve the clinically desired degree of tissue expansion. CPX™ Breast Tissue Expanders are devices that contain magnetic injection domes and are NOT MRI compatible. Do not use the CPX™ Tissue Expander in patients where an MRI may be needed. DO NOT use the CPX™ Tissue Expander in patients that have a previously implanted device that could be affected by a magnetic field. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed. Your patient needs to be informed and understand the risks and benefits of MENTOR® Tissue Expanders, and she should be provided with an opportunity to consult with you prior to deciding on surgery. For detailed indications, contraindications, warning and precautions associated with the use of all MENTOR® Implantable Devices, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at [www.Mentorwllc.eu](http://www.Mentorwllc.eu) © Mentor Worldwide LLC 2019. Intended for use by or under the direction of a physician use, it is important to read the Instructions for Use and to understand the contraindications, warnings, and precautions. All rights reserved. All other companies' brand names mentioned herein are the trademarks of their respective owners.